



# Oviduct Organoids Kit (Mouse) #BFM-001-003

## Introduction

According with the new international ethical guidelines, which suggest to reduce the number of *in vivo* experiments, new *in vitro* models that recapitulate the physiology and anatomy of living tissues are strongly requested in medical research. This kit was design to be upfront with this philosophy. 3D-organ cultures represent the future of biomedical research in which scientists could test their hypotheses in different field of interest including cancer, obesity, neuroscience, cardiovascular disease, environmental risk and nutrition.

## Safety Information

Kit contains:  
Penicillin 100 units/ml  
Streptomycin 100 µg/ml  
Gentamicin 40 µg/ml  
DMSO 0,2%

## Quality Control

Each lot of the kit was tested to ensure consistent product quality.

## Instruments

Laminar flow tissue culture hood, incubator at 37°C and 5% CO<sub>2</sub>; centrifuges for 2 and 15 ml tubes; pipettes: P1000, P200, P20 and P2; sterile tips; ice; 24 multiwell plate; automatic pipettor.

## Reagents

Name	n°	ml	storage	stability
Washing Medium (WM)	6	50	4°C	3 months
Growth Medium (GM)	10	10	-20°C	3 months
Growth Factors 1000x (GF)	10	0,01	-20°C	3 months
Recovery Sol 1 (RS1)	2	50	4°C	3 months
Recovery Sol 2 (RS2)	1	10	4°C	3 months
Optimized Extracellular Matrix (OEM)	3	0,5	-80°C	3 months
Colon Organoids (optional)	1	1	Liquid Nitrogen	1 year

**For Research Use Only (RUO).**

## **PROTOCOL**

Before starting:

- Thaw OEM overnight in ice and keep in a refrigerator. OEM is soluble in ice and solid from room temperature to 37°C.
- Add 10 µl of GF to 10 ml of GM and leave at 4°C. Do not refreeze GF and GM.

## **Thawing Organoids**

**NB:** leave GM at RT for 30' before add to organoids culture (do not pre-warm at 37°C).

- Thaw the cryovial at 37°C in a water bath by gentle agitation.
- Remove the vial from the water bath and decontaminate it by dipping in or spraying with 70% ethanol. Follow strict aseptic conditions in a laminar flow tissue culture hood for all further manipulations.
- Unscrew the cap of the vial and transfer the contents to a sterile centrifuge 15 ml tube containing 9 ml of WM.
- Spin 5' at 0.4 x g ( $\approx$  1000 rpm).
- Remove carefully the supernatant with vacuum leaving 200 µl of medium.
- Add 800 µl of WM and gently resuspend the pellet.
- Transfer in a 1.5 ml tube and centrifuge 5' at 0.4 x g
- Remove the supernatant with a p200 (do not touch the pellet).
- Resuspend the pellet in 15 µl of OEM with a sterile p20 cutted tip.
- Plate a single drop of OEM/organoids mix on the center of a well (24 mw plate).
- Invert the plate and leave in the incubator for 30'.
- Add 500 µl of GM.
- Change the GM every 2 days.
- Check the organoids every day.

## **Passage for Maintenance**

- Carefully remove the medium without touching the organoids.
- Gently wash two times with 500 µl of RS1.

- Add 200 µl of ice cold RS2 with p1000 and breakdown the OEM by pipetting up and down for 10 times.
- Incubate 5' at 37°C in incubator.
- Add 500 µl of ice cold WM and disaggregate organoids by pipetting up and down for 10 times.
- Collect the supernatant and transfer in a 15 ml tube with 9 ml of ice cold WM and leave 5' on ice.
- Centrifuge 5' at 0.4 x g.
- Remove carefully the supernatant with vacuum and leave 200 µl of medium.
- Add 800 µl of WM and gently resuspend the pellet.
- Count the organoids (only aggregate of more than 5 cells).
- Transfer in a sterile 1.5 ml tube, leave 5' on ice and centrifuge 5' at 0.4 x g.
- Plate 1000-3000 organoids in each well.

## **Ordering Information**

Code	Product	Package
BFM-001-003	Oviduct Organoids Kit (Mouse)	1 Kit

## **References**

*Sato et al.*, Single Lgr5 stem cells build crypt-villus structures in vitro without a mesenchymal niche. *Nature*. 2009 May 14;459(7244):262-5.

*Matano M et al.*, Modeling colorectal cancer using CRISPR-Cas9-mediated engineering of human intestinal organoids. *Nat Med*. 2015 Mar;21(3):256-62